



SNPase for SNP-Genotyping

- Hot-Start Polymerase for SNP detection by allele-specific PCR and micro sequencing

Features:

- 10-15 fold lower mutation rate than Taq DNA Polymerase
- high fidelity allele-specific amplification of DNA fragments
- high specificity with lowest background AS-PEX and AS-PCR
- Hot-Start activity for less primer dimers
- only 5'-3' polymerase activity, lack of 5'-exonuclease activity

Applications:

- High specific PCR
- Multiplex PCR
- Real-Time PCR with intercalation dyes
- high fidelity dNTPs and ddNTPs
- *Mini-Sequencing, SNP-genotyping*

Description:

SNPase is Taq DNA Polymerase with unique N-terminal deletion and proprietary amino acids substitutions introduced into the active center of the enzyme. This modification causes dramatic increase of sensitivity of the enzyme to mismatches at 3'-end of the primer. Consequently, non-perfect annealing of the primers does not result in unspecific amplicons formation. This enzyme has only 5'-3' polymerase activity and is recommended for SNP genotyping by allele-specific PCR (AS-PCR), allele-specific primer extension (AS-PEX) and minisequencing procedures.

Minisequencing SNP genotyping with SNPase DNA Polymerase can be carried out by the procedure described in *:

* Reference for minisequencing protocol: Lovmar L, Fredriksson M, Liljedahl U, Sigurdsson S, Syvänen AC.

Quantitative evaluation by minisequencing and microarrays reveals accurate multiplexed SNP genotyping of whole genome amplified DNA. *Nucleic Acids Res.* 2003;31:e129.

Unit definition:

One unit is defined as the amount of enzyme that incorporates 10nmoles of dNTPs into acid-insoluble form in 30 minutes at 72°C.

Concentration:

10-25 u/μl

Reaction Buffer supplied:: 5X Reaction buffer without MgCl₂,
1 tube MgCl₂ 100 mM

Note:

- optimal MgCl₂ concentration: 3.0 -3.5 mM in the 1X reaction mixture
- higher MgCl₂ concentrations results in higher yield (up to 4.5 mM)
- lower MgCl₂ (2.5 mM) results in higher specificity
- DNA fragments up to 400 bp from Human genomic DNA and 500 bp from Phage-DNA

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..



Usage:

Components	Volume per reaction
5 X reaction buffer without MgCl ₂	5 µl
MgCl ₂	2.5 - 4 mM
dNTP-Mix	0.2 mM each
primer mix (5 µM stock)	0,9-1,1 µl (5 pmol)
Template DNA	75-125 ng/25 µl genomic DNA
SNpase	0.2 - 0.5 µl (5-12 units)
Sterile dest. Water (molecular grade)	up to 25 µl total reaction volume

General Thermo-Cycler protocol:

Step	Time	Temperature
Initial denaturation	1-2 min	94-95°C
30-35 Cycles:		
Denaturation	20-30 sec	94-95°C
Annealing	15-30 sec	59-68°C
Extension	30-40 sec	68-72°C per 1kb
Final extension	5 min	72°C

Ordering Information:

Cat.-no	Description	Amount
S500	SNpase	500 units
S505	SNpase	2500 units

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..