

## Pfu/Psp DNA Polymerase

### Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'=>3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors. The enzyme has no 5'=>3' exonuclease activity.

### Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis

### Description:

Pfu/Psp DNA polymerase, isolated from the archae bacteria *Pyrococcus furiosus*/species is a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** 5 u/μl

### Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nM of dNTPs into acid insoluble material in 30 minutes at 75°C.

### Storage Buffer:

50 mM Tris-HCl, pH 8.2, 0.1 mM EDTA, 0.1% Tween 20, 0.1% Nonident P40, 1 mM DTT, 50% Glycerol

### Reaction Buffer 10 X:

100 mM KCl, 160 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20 mM MgSO<sub>4</sub>, 200 mM Tris-HCl, pH8.8, 1% Triton X-100, 1 mg/ml BSA

### Quality control:

- Tested for the DNA amplification of 2,2 kb from lambda DNA
- Contamination level check of bacterial DNA
- Purity by SDS-Page > 90 %

### Usage:

Standard protocol:

- Do not use dUTP or dITP or primers containing these nucleotides

Components	Volume per reaction	end conc.
10X reaction buffer with MgSO <sub>4</sub>	5 μl	1X
dNTP-Mix (40mM = 10mM each)	1.0 μl	200 μM each
Up-stream primer (e.g. 20 μM)	0,5 μl	0.1-1.0 μM
Down-stream primer (e.g. 20 μM)	0.5 μl	0.1-1.0 μM
Template DNA (10 ng/μl)	1.0 μl	<= 0,5 μg
Pfu/Psp DNA Polymerase (5 u/μl)	0.2 - 0,4 μl	1-2 units
Sterile dest. Water (molecular grade)	up to 50 μl	

### Note:

- vortex all solutions carefully before using

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

*.. a good decision ..*

dispense all reagents on ice to avoid degradation of primers and dNTP's  
 - add the enzyme after Template DNA  
 - may you have to optimize the MgSO<sub>4</sub> concentration for best result

**General Thermo-Cycler protocol:**

Step	Time	Temperature
Initial denaturation	1-3 min	95°C
<b>25-35 Cycles:</b> Denaturation Annealing Extension	30-100 sec 30-65 sec 1-2 min (per 1kb)	95°C 37-69°C 72-75°C
Final extension	5 min	72-75°C

**Storage:** at -20°C for 24 months

**Transportation:** on blue ice

**Related products:**

**Ordering information:**

Cat.-no	Description	Amount
S116	Pfu/Psp DNA Polymerase	1x250 units
S117	Pfu/Psp DNA Polymerase	2x250 units
S118	Pfu/Psp DNA Polymerase	10x250 units

*.. a good decision ..*