

GelRED™

Nucleic Acid Gel Stain – 10000X

Product description

GelRed™ is an ultra-sensitive, extremely stable and environmentally safe fluorescent nucleic acid dye designed to replace the highly toxic ethidium bromide (EB) for staining dsDNA, ssDNA or RNA in agarose gels or polyacrylamide gels. GelRed is far more sensitive than EB without requiring a destaining step. GelRed and EB have virtually the same spectra, so you can directly replace EB with GelRed without changing your existing imaging system.

GelRed™ can be used to stain dsDNA, ssDNA or RNA in agarose gel via either precast or post gel staining. GelRed can also be used to stain dsDNA, ssDNA or RNA in polyacrylamide gel via post gel staining. GelRed is also compatible with downstream DNA manipulations such as digestion with a restriction enzyme, Southern blotting techniques and cloning.

A series of safety tests have confirmed that GelRed™ is noncytotoxic, nonmutagenic and nonhazardous at concentrations well above the working concentrations used in gel staining. As a result, GelRed can be safely disposed of in regular trash, providing convenience and reducing cost in waste disposal. GelRed™ is supplied as 10,000X solution in water or for your convenience.

As nucleic acid binding dyes can affect DNA migration during electrophoresis, post-staining of gels is highly recommended. Post-staining with GelRed™ results in superior sensitivity and eliminates the possibility of dye interference with DNA migration. Post-staining with GelRed™ is simple, requiring no destaining and no special buffer. Simply dilute the concentrated dye in 0.1M NaCl or water and incubate the gel in the diluted dye solution for 30 minutes. The staining solution is perfectly stable at room temperature, permitting it to be used multiple times, and is substantially more sensitive than that using EB. Although the post-staining method is recommended, precast gels may also be tried with GelRed™. However, some DNA samples, such as those derived from plasmid DNA digestion by certain

restriction enzymes, may experience migration retardation or compromised resolution. Thus, both the post-stained and precast gels can be performed to determine which one may better meet your needs. GelRed™ can also be used to stain dsDNA*, ssDNA or RNA in polyacrylamide gel via post gel staining. Precast polyacrylamide gel staining with GelRed™ is not recommended because of relatively high background fluorescence.

* GelRed™ is twice as sensitive to dsDNA than ssDNA or RNA.

Gel staining with GelRed™ is compatible with downstream DNA manipulation such as digestion with a restriction enzyme, Southern blotting techniques and cloning. GelRed™ may be removed from DNA by ethanol precipitation.

GelRed™ Nucleic Acid Gel Stain – 10000X in water is a concentrated solution that can be diluted 10,000 times for use in precast gel staining or 3,300 times for use in post gel staining according to the procedures described.

500µL of 10000X solution can be used to stain at least 100 precast gels or post-stain at least 100 minigels

Staining Protocols

1. Staining DNA by Post Gel Staining

a. Run gels as usual according to your standard protocol.

b. Dilute the GelRed™ 10000X stock solution about 3,300 fold to prepare a 3X staining solution in water with 0.1M NaCl (e.g. add 15µL of GelRed™ stock solution and 5mL NaCl to 45mL water).

While GelRed™ 1X staining solution can also be used for post gel staining, the sensitivity is generally less than with 3X staining solution

(NOTE: inspect the 10000X vial carefully. If dye precipitation occurs, heat up or sonicate the vial).

Use of NaCl in the staining solution is optional. NaCl in the solution enhances the staining, but may promote dye precipitation if the staining solution is to be used repeatedly. Any staining solution to be re-used is preferably stored at room temperature in a dark place to reduce possible dye precipitation problems.

c. Carefully place the gel in a suitable container such as a polypropylene container. Gently add sufficient amount of the 3X staining solution to submerge the gel.

d. Agitate gel gently at room temperature for about 30 minutes. Optimal staining time may vary somewhat depending on the thickness of the gel and the percentage of agarose. For polyacrylamide gels containing 3.5 – 10% acrylamide, typical staining time is 30 minutes up to 1 hour with gels of higher acrylamide content requiring longer staining time. The staining solution can be re-used at least 2-3 times.

The unused staining solution can be stored at room temperature in a dark place.

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462
E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..

e. View the stained gel with a standard transilluminator (302nm or 312nm) and photograph the gel using an ethidium bromide filter. Similarly, a SYBR(R) or GelStar(R) filter may also be used for photographing with equally good results.

2. Staining of DNA by Precasting GelRed™ Gels*

a. Prepare agarose gel solution using your standard protocol.

b. Dilute GelRed™ 10000X stock solution into the agarose gel solution at 1:10000 (e.g. 5µL stock solution added to 50mL of gel solution). Since GelRed™ is generally thermally stable the 10000X stock solution can be added while the gel solution is still hot. Make sure that the dye is thoroughly mixed with the gel solution by swirling, stirring, or inversion (**NOTE:** inspect the 10000X vial carefully. If dye precipitation occurs, heat up or sonicate the vial).

Alternatively the stock solution may be pre-combined with agarose powder and electrophoresis buffer of your choice followed by microwaving or other heating procedures commonly used for preparing agarose gels. GelRed™ is compatible with all commonly used electrophoresis buffers.

c. Cast the gel and allow it to solidify. Any leftover gel solution may be stored and re-heated later for additional gel casting. Since GelRed™ is hydrolytically stable, GelRed™ precast gels may be prepared in large quantities and stored for later use. To avoid mould formation, we recommend that the precast gels be stored in a refrigerator.

d. Load samples and run the gels using your standard protocol.

NOTE: Use only 1/3 or less of your normal DNA-marker amount to avoid overexposure (no bands of marker visibly separated).

e. View the stained gel with a standard transilluminator (302nm or 312nm) and photograph the gel using an ethidium bromide filter. Similarly, a SYBR(R) or GelStar(R) filter may also be used for photographing with equally good results. **NOTE:** If you consistently see band smearing and/or poor band separation, run a gel and post-stain by following the protocol provided (point 1) to confirm if the problem is caused by the dye or other factors unrelated to the dye. If post gel staining is normal and the problem is not caused by the dye, try any of the followings:

- Lower the amount of nucleic acid loaded.
- Lower running voltage
- Lower the amount of agarose in the gel
- Run a longer gel
- Increase the thickness of the gel
- Increase the solidification time to ensure sharp well formation
- Improve your loading technique or select post gel staining

* **Precasting GelRed™ gel is not suitable for acrylamide gels. Use post gel staining for acrylamide gels.**

Features of GelRed™

Safer than EB

Shown by the Ames test and other tests to be nonmutagenic and noncytotoxic

Easy disposal

Passed environmental safety tests for direct disposal down the drain or in regular trash

Ultra-sensitive

Much more sensitive than EtBr and SYBR Safe

Extremely stable

Available in water, stable at room temperature for long-term storage and microwavable

Simple to use

Very simple procedures for either precast and post gel staining

Perfectly compatible with a standard UV transilluminator

GelRed replaces EtBr with no optical setting change

Perfectly compatible with downstream applications

Compatible with downstream DNA manipulations such as digestion with a restriction enzyme, Southern blotting techniques and cloning.

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..




Figure 1. GelRed™ is significantly more sensitive than ethidium bromide (EB) for detecting low-level DNA, especially in the lower molecular weight area. Shown left are two-fold serial dilutions of 1 Kb Plus DNA Ladder from Invitrogen electrophoresed on 1% agarose gels precasted with GelRed or EB in 1x TBE. The total amount of DNA loaded per lane was: 200ng, 100ng, 50ng and 25ng from left to right. Gels were imaged using 300 nm transillumination and photographed with an EB filter and Polaroid 667 black-and-white print films.




Figure 2. Excitation and emission spectra of GelRed and GelGreen in the presence of DNA in PBS buffer

Storage and Stability

GelRed™ is a very stable dye. We recommend that the 10000X solutions are stored at RT in a dark place. The solution may also be stored at a lower temperature such as 4°C. Dye precipitation may occur during prolonged low temperature storage. When this occurs, heat up the solution in a hot water bath at 45°C to 50°C for two minutes and vortex the solution.

The 1X or 3X post column staining solutions may also be stored at room temperature in a dark place for at least one year. Exposure to light should be avoided during long-term storage. However, the dye can be handled under ambient light without any problem during staining experiments.

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..




Figure 3. Normalized absorbances of GelRed™ and SYBR® Gold 1X TBE gel-staining solutions at 500 and 488nm respectively over time at room temperature. The starting absorbance values for GelRed™ and SYBR® gold were 0.029 and 0.051, respectively.

Note: *GelRed and its uses are covered by pending US and international patents.
 **SYBR is trademark of Molecular Probes, Inc.

Toxicity

GelRed™ was subjected to a series of tests by three independent testing services to assess the dye's safety for routine handling and disposal.

These tests include:

1. Glove penetration test.
2. Cell membrane permeability and cytotoxicity test.
3. Ames test.
4. Environmental safety tests. Test results confirm that the dye is impenetrable to both latex gloves and cell membranes. The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. GelRed™ appears to be completely cell membrane-impermeable, which may be a key factor responsible for the observed low toxicity. However, since these tests were not performed on human, we still advise that researchers exercise precautions when handling the dye or any other DNA-binding molecules by wearing protective gears.

Disposal

GelRed™ has successfully passed environmental safety tests in compliance with CCR Title 22 Hazardous Waste Characterization. As a result, GelRed™ is not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash, providing convenience and reducing cost in waste disposal.

Note: *GelRed and its uses are covered by pending US and international patents.
 **SYBR is trademark of Molecular Probes, Inc.

Trouble shooting

Problem Suggestion

- Smeared DNA bands in gel 1.Reduce the amount of DNA loaded by one-half to one-third.
 Blown out or smeared bands can be caused by overloading. This is frequently observed with DNA ladders.
- 2.Perform post-staining instead of pre-casting.
 - 3.Pour a lower percentage agarose gel for better resolution or large fragments.
 - 4.Change the running buffer. TBE buffer has a higher buffering capacity than TAE
- Discrepant DNA migration in pre-cast gel
 GelRed is designed to be larger than other dyes to prevent it from entering cells, thus rendering the dye safer. The migration of DNA may be affected depending on the dye:DNA ratio.
- 1.Reduce the amount of DNA loaded by on-half to one-third.
 - 2.Reduce the amount of dye used, i.e. use 0.5X in pre-cast gels.
 - 3.Post-stain the gel in 3X GelRed to avoid any interference the dye may have on migration during electrophoresis.

.. a good decision ..

Weak fluorescence, decreased dye performance over time, or film of dye remains on gel after post-staining. The dye may have precipitated out of solution.

- 1-Heat GelRed solution to 45-50°C for two minutes and vortex to redissolve.
- 2.Store dye at room temperature to avoid precipitation.

Frequently Asked Questions

Question Answer

Can GelRed be used to stain ssDNA or RNA?

GelRed can be used to stain ssDNA and RNA. GelRed is 5-times more sensitive for single stranded nucleic acids than GelGreen.

Is GelRed compatible with downstream applications such as cloning, ligation and sequencing?

Yes. We recommend gel extraction kits, Exo-Sap protocol or phenol-chloroform extraction to remove the dye from DNA.

Is GelRed compatible with Southern or northern blotting?

Yes. GelRed can be used for blotting.

Can I reuse GelRed pre-cast gel after electrophoresis?

Yes, it is possible, but we do not recommend reusing GelRed pre-cast gels as signal decreases with subsequent electrophoresis.

How should I dispose GelRed?

GelRed has passed EPA regulated Title 22 test. Some facilities have approved the disposal of GelRed directly down the drain. However, because regulations vary, please contact your safety office for local disposal guidelines. GelRed can be adsorbed to activated carbon (also known as activated charcoal) for disposal as chemicals waste.

What is the lower detection limit of GelRed?

Some users have reported being able to detect less than 0.1ng DNA. However, the limit of detection will depend on instrument capability and exposure settings.

Does GelRed need to be used in the dark?

You can use the dye in room light, however we recommend storing the dye in the dark.

Is there a difference between 10,000X GelRed in DMSO and water?

The GelRed stock in water is a newer and improved product compared to the stock in DMSO. We recommend using GelRed in water to avoid the potential hazards of handling DMSO, a solvent that can be absorbed through the skin.

Ordering information:

Cat.-No:	Description	Amount
S420	GelRED™	500 µl
S425	GelRED™	5x500µl

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..