

Specification

Product	LB Medium powder according to Lennox for molecular biology
Brand	GeneOn GmbH
Product number	CE182, ff.
Synonym(s)	Luria-Bertani Medium, low salt
Storage	below 25°C
Retest period (months)	48
For cultivation and maintenance of recombinant strains of Escherichia coli and other not particularly fastidious microorganisms.	
Parameter	Specification
Appearance	cream to yellow coloured free flowing powder
pH value	6.80 - 7.20
Composition:	
Yeast extract	5 g/l
NaCl	5 g/l
Tryptone (Casein enzymic hydrolysate)	10 g/l
Cultural response (E. coli)	luxuriant