

Product	BCIP for molecular biology
Product number	CE108 ff.
Synonym(s)	5-Bromo-4-chloro-3-indolyl phosphate p-Toluidine salt, X-Phos p-Toluidine salt, X-Phos-pT
Formula	$C_8H_6BrClNO_4P \cdot C_7H_9N$
Molecular weight	433.64 g/mol
CAS number	6578-06-9
Storage	-20°C
Retest period (months)	48
Parameter	Specification
Appearance	white to light brown powder
Assay	min. 98 % (revised)
DNases/RNases/Proteases	not detectable
Identity	passes test
Water (K.F.)	max. 3 % (revised)
Ludwigshafen	16 May 2023