

<b>Product</b>	<b>NADP sodium salt for molecular biology</b>
Brand	GeneON-BioScience
Product number	CE200 ff.
Synonym(s)	β-Nicotinamide adenine dinucleotide phosphate monosodium salt
Formula	C <sub>21</sub> H <sub>27</sub> N <sub>7</sub> NaO <sub>17</sub> P <sub>3</sub>
Molecular weight	765.39 g/mol
CAS number	1184-16-3
Storage	-20°C
Retest period (months)	36
<b>Parameter</b>	<b>Specification</b>
DNases/RNases/Proteases	not detectable
Assay (enzym.)	min. 95 %
Water (K.F.)	max. 6.0 %
UV Absorption	
A250/A260 (pH 7.5)	0.83 ± 0.03
A280/A260 (pH 7.5)	0.21 ± 0.02
Ludwigshafen	5 June 2022