

Specification

Product	Formamide deionized for molecular biology
Product number	CE145, ff.
Synonym(s)	Carbamaldehyde, Formic acid amide, Methanamide
Formula	CH ₃ NO
Molecular weight	45.04 g/mol
CAS number	75-12-7
Storage	2 - 8°C
Retest period (months)	48
Parameter	Specification
Appearance	clear, colorless liquid
DNases/RNases/Proteases	not detectable
Assay (from N)	min. 99.5 %
Water (non-deionized)	max. 0.1 % (revised)
Chloride	max. 0.00005 %
Fe	max. 0.00001 %
Pb	max. 0.00001 %